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Method of Test for Pesticide Residues in Foods-Multiresidue Analysis (5)

Taiwan Food and Drug Administration

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1. Scope: The test method is suitable for 310 pesticide Multi-residues analysis in food from fruits and vegetables, cereals, dried beans, tea, spice plants and other herbs etc.

2. Test methods: the sample is pre-treated using QuEChERS method(Quick-Easy-Cheap-Effective-Rugged-Safe) and then for the method of analysis with Liquid Chromatography Tandem Mass Spectrometry LC/MS/MS) and Gas Chromatograph/Tandem mass spectrometer, GC/MS/MS).

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2.1.1 Liquid Chromatograph/Tandem Mass Spectrometer:

2.1.1.1 Ion Source: Electrospray Ionization, ESI).

2.1.1.2 Chromatography Column: Acquity UPLCR HSS T3, 1.8um, an inner diameter of 2.1 mm × 10 cm, or products with same quality.

2.1.2 Gas Chromatography Tandem Mass Spectrometry:

2.1.2.1 Ion Source: Electron Impact Ionization, EI.

2.1.2.2 Chromatography Column: DB-5MS UI capillary, inner film layer thickness 0.25um, 0.25 mm ID × 30 m, or products with same quality..

2.1.3. Stirring homogenizer (Blender).

2.1.4 Grinder.

2.1.5 High-speed oscillation grinding tissue homogenizer (SPEX SamplePrep 2010 GenoGrinderR): 1000 rpm or more, or products with same quality.

2.1.6 Centrifuger: Up to 3000 × g or more, control the temperature up to those below 15°C.

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2.2 Reagent:

Acetic acid,

Formic acid and ammonium acetate are of GR/Guaranteed Reagent Grade;

Hexane and Acetone are of residual grade;

Acetonitrile and methanol are of Liquid Chromatography Grade.

Anhydrous sodium acetate, Anhydrous magnesium sulfate, Primary secondary amine(PSA),

Octadecylsilane, end-capped(C18 EC) and Graphitized Carbon Black(GCB) are of analytical grade;

Deionized water (resistivity at 25°C Up to 18 M ohm · cm or above);

310 pesticide Reference Standards as 3-keto carbofuran (items shown in Table I. Table 2 and Table 3);

Triphenyl phosphate (triphenylphosphate, TPP) internal standard.

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2.3.2 Membrane: Pore size 0.22um, PVDF material.

2.3.3 Volumetric flasks: 25 mL and 50 mL, amber.

2.3.4 Homogeneous ceramic stone (Ceramic homogenizer) (Note 1): Using Bond Elut QuEChERS P/N 5982-9313, or products with same quality.

2.3.5 Extracting powder (Note 2): with Anhydrous magnesium sulfate ($MgSO_4$) 4 g, and Anhydrous sodium acetate (NaOAc) 1 g.

2.3.6 Purification centrifuge tube Type I (Note 2): with PSA 300 mg and Anhydrous magnesium sulfate ($MgSO_4$) 900 mg, extracting solution loading capacity 6 mL, suitable for fruits and vegetables samples with high water content.

2.3.7 Purification centrifuge tube Type II (Note 2): with PSA 300 mg, C18EC 300 mg and Anhydrous magnesium sulfate ($MgSO_4$) 900 mg, extracting solution loading capacity 6 mL, suitable for cereals samples with high wax, fat and sugar content.

2.3.8 Purifying centrifuge tube Type III (Note 2): with PSA 450 mg, Anhydrous magnesium sulfate ($MgSO_4$) 900 mg, C18 EC 300 mg and GCB 50 mg, extracting solution loading capacity 6 mL, suitable for samples with high pigment content and tea analogues.

Note 1: The Ceramic homogenizer could be arbitrarily used according to the user's assessment of the viscosity of sample.

Note 2: According to their demands, the users could make their assessment for applying the various extraction and purification kits commercially available.

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2.4 Mobile Phase Solution Preparations:

2.4.1 Mobile Phase Solution A:

Take 50 mL of Methanol and mix with 450 mL of deionized water, add and to dissolve Ammonium Acetate 0.19g

and mixed completely, filter with membrane. The filtrate is for mobile phase A.

2.4.2 Mobile Phase Solution B:

Take 450 mL of Methanol and mix with 50 mL of deionized water, add and to dissolve Ammonium Acetate 0.19g

and mixed completely, filter with membrane. The filtrate is for mobile phase B.

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2.5 Internal Standard solution preparation:

Take Triphenylphosphate(TPP) internal standard about 40 mg, accurately weighed, dissolved in methanol and dilute to 50 mL, as an Internal Standard Stock Solution, stored at -18°C in dark for use.

2.5.1 Take appropriate amount of Internal Standard Stock Solution and then dilute with Methanol to 75 μ g/mL, that will be the Internal Standard Solution used for the Test Solution Preparation in section 2.8.

2.5.2 Take appropriate amount of Internal Standard Stock Solution and then dilute with Methanol to 7.5 μ g/mL, that will be the Internal Standard Solution used for LC/MS/MS analysis in section 2.9.1

2.5.3 Take appropriate amount of Internal Standard Stock Solution and then dilute with Acetone to 7.5 μ g/mL, that will be the Internal Standard Solution used for GC/MS/MS analysis in section 2.9.2

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2.6 Reagents Preparation:

2.6.1 1% Acetic Acid(HOAc) in Acetonitrile(MeCN) solution:

Take 10 mL of Glacial Acetic Acid(HOAc) and mixed with 990 mL of Acetonitrile(MeCN).

2.6.2 5% Formic Acid in Acetonitrile(MeCN) solution:

Take 5 mL Formic acid and mixed with 95 mL of Acetonitrile(MeCN).

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2.7 Standard Solution Preparation:

2.7.1 Take about 25 mg pesticide reference standards, accurately weighed and dissolved in Acetonitrile(MeCN) to 25 mL

as the Standard Stock solution, store at -18 $^{\circ}$ C in dark for use.

Take appropriate amount of the Standard Stock solution and diluted with Methanol to 1 μ g/mL, as the Standard solution used for LC/MS/MS analysis in section 2.9.1.

2.7.2 Take about 25 mg pesticide reference standards, accurately weighed and dissolved in Acetone or Hexane to 25 mL

as the Standard Stock solution, store at -18 $^{\circ}$ C in dark for use.

Take appropriate amount of the Standard Stock solution and diluted with Acetone: Hexane (1:1, v/v) solution

to 1 µg/mL,

as the Standard solution used for GC/MS/MS analysis in section 2.9.2.

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2.8.1 Fruits and vegetables, spice plants and other herbs(fresh meal):

Accurately weigh about 10g homogeneous sample into a centrifuge tube, frozen and then add 1% Acetic Acid(HOAc) in Acetonitrile(MeCN) solution 10 mL, and 75µg/mL Internal Standard Solution 10µL, and then sequentially add one Ceramic homogenizer and extracting powder.

Close the cap and immediately followed with intense oscillations for times to prevent from caking of the salt, and then vigorously shaken by hand or by high speed tissue grinding and oscillating homogenizer in 1000rpm for 1 minute, and centrifuged at 3000 × g for 1 min at 15°C.

Transfer the supernatant 6 mL to a Purification centrifuge tube(TypeI) and then vigorously shaken by hand or by high speed tissue grinding and oscillating homogenizer in 1000rpm for 1 minute, and centrifuged at 3000 × g for 2 min at 15°C.

Transfer the supernatant 1 mL to a scaled vials and evaporated with nitrogen to just dried,

and the residue was dissolved in an appropriate amount of Methanol, with 5% Formic Acid in Acetonitrile(MeCN) solution 10 μ L added, to get 1mL volume. Mix homogenously and filter with membrane to get the Test Solution I, for the LC/MS/MS analysis.

Transfer another supernatant 1 mL to a scaled vials and evaporated with nitrogen to just dried, and the residue was dissolved in an appropriate amount of Acetone: Hexane (1: 1, v/v) solution, with 5% Formic Acid in Acetonitrile(MeCN) solution 10 μ L added, to get 1mL volume. Mix homogenously and filter with membrane to get the Test Solution II, for the GC/MS/MS analysis.

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2.8.2 Cereals and dried beans:

Accurately weigh about 5g homogeneous grinded sample into a centrifuge tube, add pre-cooling deionized water 10 mL and stand for 20 minutes, and then add 1% Acetic Acid(HOAc) in Acetonitrile(MeCN) solution 10 mL. The following procedures are the sample as that in the section 2.8.1, except that the Purification centrifuge tube(Typell) is used instead.

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2.8.3 Tea, spice plants and other herbs (dried):

Accurately weigh about 2g homogeneous grinded sample into a centrifuge tube,

add pre-cooling deionized water 10 mL and stand for 20 minutes,

and then add 1% Acetic Acid(HOAc) in Acetonitrile(MeCN) solution 10 mL.

The following procedures are the same as that in the section 2.8.1,

except that the Purification centrifuge tube(Typelll) is used instead.

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2.9 Matrix-matched calibration curve preparation

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2.9.1 LC/MS/MS:

Take a blank sample and process according to the 2.8.Test Solution Preparation steps without the addition of Internal Standard solution to produce the supernatant after SPE purification.

For each standard levels, transfer 1mL supernatant to scaled vial and evaporate with nitrogen to just dried,

and the residue was dissolved in an appropriate amount of Methanol, with 1μg/mL (Note 3) Standard Solution 2 ~ 200μL respectively,

7.5μg/mL of Internal Standard solution 10μL, and 5% Formic Acid in Acetonitrile(MeCN) solution 10μL added,

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The Matrix-matched calibration curves are produced as a plot of (peak area/internal standard peak area)ratio of each pesticides versus the corresponding concentration of each pesticides in the 0.002 ~ 0.2 μ g/mL range (Fipronil is 0.0004 ~ 0.04 μ g/mL).

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Liquid Chromatography Tandem Mass Spectrometry measurement conditions (Note 4):

Chromatography Column: Acquity UPLCR HSS T3, 1.8 μ m, 2.1mm(ID) \times 10 cm.

Mobile phase solution: Analysis with the following Gradient conditions for A and B solution

Mobile phase flow rate: 0.3 mL/min.

Capillary voltage :

Positive ion electrospray ionization (ESI +) using 3.5 kV,

Negative ion electrospray ionization (ESI-) using 1.6 kV.

Ion source temperature (Ion source temperature): 150 $^{\circ}$ C.

Solvent volatilization temperature (Desolvation temperature): 450 $^{\circ}$ C.

Injection cone gas flow rate (Cone gas flow): 30 L/hr.

Flow rate of vaporized solvent (Desolvation flow): 900 L/hr.

Detection mode: Multiple reaction monitoring(MRM).

Detective Measured ion pair, injection cone voltage (cone voltage) and collision energy (collision energy) as shown in Table 1 and Table 2.

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2.9.2 GC/MS/MS:

Take a blank sample and process according to the 2.8.Test Solution Preparation steps without the addition of Internal Standard solution to produce the supernatant after SPE purification. For each standard levels, transfer 1mL supernatant to scaled vial and evaporate with nitrogen to just dried, and the residue was dissolved in an appropriate amount of Acetone: Hexane(1:1, v/v)solution, with $1\mu\text{g/mL}$ Standard Solution 4 ~ 500 μL respectively, $7.5\mu\text{g/mL}$ of Internal Standard solution 10 μL , and 5% Formic Acid in Acetonitrile(MeCN) solution 10 μL added, to get 1mL volume. mix homogenously.

The prepared Matrix-matched calibration standards are analyzed by the following conditions.

The Matrix-matched calibration curves are produced as a plot of (peak area/internal standard peak area)ratio of each pesticides

versus the corresponding concentration of each pesticides in the 0.004~0.5 $\mu\text{g/mL}$ range.

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Gas Chromatography Tandem Mass Spectrometry measurement conditions (Note 4):

Chromatography Column: DB-5MS UI capillary, inner film layer thickness 0.25um, 0.25 mm ID × 30 m.

Chromatography Column temperature: initial temperature: 60°C, 1 min;

Heating rate: 40°C/min;

Temperature: 170°C;

Heating rate: 10°C/min;

Final temperature: 310°C, 2.25 min.

Mobile phase flow rate: Helium, 1 mL/min.

Injector temperature : 280°C.

Injection mode: splitless.

Ionization modes: electron impact, 70 eV.

Ion source temperature: 300°C.

Detection mode: Multiple reaction monitoring, MRM, as detection ion pairs and collision energy shown in Table III.

Note 3: Select appropriate standard solution spikings for producing the Matrix-matched calibration curve of Fipronil.

Note 4: Setting suitable measurement conditions according to the instrument used if the measurement conditions listed above are not suitable.

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2.10 Identification Test and Quantitation:

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2.10.1 LC/MS/MS: ,

Exactly measure 10μL of Test solution and Standard solution and inject into the Liquid Chromatography Tandem Mass spectrometry respectively.

Perform analysis in accordance with the measurement conditions from section 2.9.1,

Identifications are based on the peak retention time and the Multiple reaction monitoring(MRM) relative ion

strength

of the Test solution and Standard solution(Note 5), and calculated in accordance with the following formula, to obtain the content of each pesticides(ppm) of samples :

The content of each pesticide in a sample (ppm) = $C \times V / M$

C: Each pesticides concentrations($\mu\text{g}/\text{mL}$) of Test solution obtained/indexed from corresponding Matrix-Matched calibration curve for each pesticides.

V: Volume of extracted sample(10 mL) with 1% Acetic Acid(HOAc) in Acetonitrile(MeCN) solution.

M: Sampled weight of the sample for analysis(g)

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2.10.2 GC/MS/MS:

Exactly measure $10 \mu\text{L}$ of Test solution and Standard solution and inject into the Gas Chromatography Tandem Mass spectrometry respectively.

Perform analysis in accordance with the measurement conditions from section 2.9.2,

Identifications are based on the peak retention time and the Multiple reaction monitoring(MRM) relative ion strength

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Note 5: Relative Ionic Strength is obtained from the from the ratio of the peak areas from the qualitative and quantitative ion pair ($\leq 100\%$),

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1. The Quantitation Limits of this test method are shown in Table I, Table 2 and Table 3.
2. The items listed in this test method could be assessed for the demands to the analysis with GC/MS/MS or LC/MS/MS.
3. The method does not apply to analysis of Benfuracarb and Pymetrozine in Tea matrix.
4. If there are interfering substances in the food that affect the result of analysis, the practioner should check and study the cause and solution.

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TableI, The Multiple reaction monitoring(MRM) parameters and Quantitative detection limits(by LC/MS/MS in positive ion mode) for 144 pesticide items as 3-keto carbofuran and Internal Standard.

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TableII, The Multiple reaction monitoring(MRM) parameters and Quantitative detection limits(by LC/MS/MS in negative ion mode) for 6 pesticide items as Talon.

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TableIII, The Multiple reaction monitoring(MRM) parameters and Quantitative detection limits(by LC/MS/MS in positive ion mode) for 160 pesticide items as Acetochlor and Internal Standard.

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a Suitable for vegetables and fruit, spice plants and other herbs(fresh meal).

b Applies to cereals and dried beans.

c Suitable for tea, spice plants and other herbs(dried).