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GL-LF096 Fraction Collector adopts the latest ultra-high speed 32-bit RISC central processor computing and 32-bit DSP processor for core control.

It is based on the all-digital frequency generation system, advanced electromagnetic valve components and friendly user interface, making the operation and control of the device very convenient.

GL-LF096 fraction collector using independent valve module control · collection method and collecting channel can be set according to peak height, peak time.

Also, can also be controlled via the front panel control keys for channel control 。A variety of communication modes selection, as standard RS232, RS485, Ethernet interface · are available for GL-LF096 fraction collector control.

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
The Power supply unit is composed of control power supply and driving power supply sub-unit, which provides 24V control power for the GL-LF096 fraction collector.

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CPU controller unit is a core part of the fraction collector, which allows the user to select the fraction collection working methods, to edit time and other programs.

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XY module is a fluid control device, it ensures that the liquid fraction can smoothly and accurately reach the entrance of each tube , enabling smoothly preparation works.

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The PC monitor is used for Chromatogram Workstation software control, and the operation is simple and intuitive 。 Also, the 7" true color display is used as spare tool for program editing and control if necessary.

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Unpacking the shipped instrument and check if any damage; if there is any damage during transport, ask statement as soon as possible to the instrument carrier sector for the relevant payment. According to accompanying packing list for the shipping, check to see

if all corresponding equipment and accessories are included.

If the shortage, please contact your dealer or by directing a request to the manufacturers. (See packing list)

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Check the working room for GL-LF096 fraction collector, and the environment should meet :

供電電壓 Supply voltage	115V-230VAC · 接地電阻 < 2Ω 115V-230VAC, ground resistance < 2Ω
環境溫度 Ambient temperature	常溫 Room temperature
環境濕度 Environment humidity	20%-80%
安裝場地 Installation site	安裝場所應遠離強震動源及強磁場干擾 Installation site should be away from strong vibration source and strong magnetic field

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According to the user's needs for different preparations，install FEP tubes that are enough by passing 200mL solvent per minute.

Tubes installation methods:

Remove the fastening screws on the fraction mouth, then insert tube head into a ferrule(tightening ring part) and connect the tube head firmly into position.

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After completing 3.1 and 3.2 of the steps, insert to connect the power cord into the instrument rear panel power input socket, and connect to power source.

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The rear panel diagram of GL-LF096 fraction collector.

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Make sure that the instrument power cord is plugged to a power source socket with 3 pins that is in correct voltage and loading capacity, and is well grounded.

Wrong power supply may cause serious personal injury.

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When first running the fraction collector, please run in small flow rate with careful monitoring. If there is no liquid carried within few minutes of continuous operation, shut down immediately to check whether the pump is faulty, and it must be resolved before to power-on again.

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Turn on the power, the fraction collector boots into Home Page.

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


## Second Boot interface of the GL-LF096 fraction collector

## 4.2 工作站軟體操作 Workstation software operations

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Click Copy button to get fraction data.

Fraction

Status												
Motion Speed	Collect Flowrate	Tube Volume (%)	Start Time (Min)	Start Time (Sec)	Operation Time Interval (Min)	Operation Time Interval (Sec)	Cycle Time	Zero Tube Diameter	Zero Tube Length	One Cycle Tube Number	End Time (Min)	End Time (Sec)
32.1	12.34	10.5	1	30	2	30	1	5	9	26	4	0

Setting												
Motion Speed	Collect Flowrate	Tube Volume (%)	Start Time (Min)	Start Time (Sec)	Operation Time Interval (Min)	Operation Time Interval (Sec)	Cycle Time	Zero Tube Diameter	Zero Tube Length	One Cycle Tube Number	End Time (Min)	End Time (Sec)
0	0	0	0	0	0	0	0	0	0	0	0	0

Run / Stop												
<input type="button" value="Run"/>												<input type="button" value="Copy"/> <input type="button" value="Apply"/>
												<input type="button" value="Close"/>

Fraction

Status												
Motion Speed	Collect Flowrate	Tube Volume (%)	Start Time (Min)	Start Time (Sec)	Operation Time Interval (Min)	Operation Time Interval (Sec)	Cycle Time	Zero Tube Diameter	Zero Tube Length	One Cycle Tube Number	End Time (Min)	End Time (Sec)
32.1	12.34	10.5	1	30	2	30	1	5	9	26	4	0

Setting												
Motion Speed	Collect Flowrate	Tube Volume (%)	Start Time (Min)	Start Time (Sec)	Operation Time Interval (Min)	Operation Time Interval (Sec)	Cycle Time	Zero Tube Diameter	Zero Tube Length	One Cycle Tube Number	End Time (Min)	End Time (Sec)
32.1	12.34	10.5	1	30	2	30	1	5	9	26	4	0

Run / Stop												
<input type="button" value="Run"/>												<input type="button" value="Copy"/> <input type="button" value="Apply"/>
												<input type="button" value="Close"/>

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Set parameters value and click Apply button.

Fraction

Status												
Motion Speed	Collect Flowrate	Tube Volume (%)	Start Time (Min)	Start Time (Sec)	Operation Time Interval (Min)	Operation Time Interval (Sec)	Cycle Time	Zero Tube Diameter	Zero Tube Length	One Cycle Tube Number	End Time (Min)	End Time (Sec)
32.1	12.34	10.5	1	30	2	30	1	5	9	26	4	0

Setting												
Motion Speed	Collect Flowrate	Tube Volume (%)	Start Time (Min)	Start Time (Sec)	Operation Time Interval (Min)	Operation Time Interval (Sec)	Cycle Time	Zero Tube Diameter	Zero Tube Length	One Cycle Tube Number	End Time (Min)	End Time (Sec)
32.1	12.34	10.5	2	30	3	30	1	5	9	26	6	0

Run / Stop												
<input type="button" value="Run"/>												<input type="button" value="Copy"/> <input type="button" value="Apply"/>
												<input type="button" value="Close"/>

#### 4.2.3.3.完成編輯 Done editing

Fraction													
Status													
Motion Speed	Collect Flowrate	Tube Volume (%)	Start Time (Min)	Start Time (Sec)	Operation Time Interval (Min)	Operation Time Interval (Sec)	Cycle Time	Zero Tube Diameter	Zero Tube Length	One Cycle Tube Number	End Time (Min)	End Time (Sec)	
32.1	12.34	10.5	2	30	3	30	1	5	9	26	6	0	

Setting													
Motion Speed	Collect Flowrate	Tube Volume (%)	Start Time (Min)	Start Time (Sec)	Operation Time Interval (Min)	Operation Time Interval (Sec)	Cycle Time	Zero Tube Diameter	Zero Tube Length	One Cycle Tube Number	End Time (Min)	End Time (Sec)	
32.1	12.34	10.5	2	30	3	30	1	5	9	26	6	0	

Run / Stop													
<input type="button" value="Run"/>													<input type="button" value="Close"/>
<input type="button" value="Copy"/> <input type="button" value="Apply"/>													

#### 4.2.4 運行餾分: Run Fraction

Click Run button

Fraction													
Status													
Motion Speed	Collect Flowrate	Tube Volume (%)	Start Time (Min)	Start Time (Sec)	Operation Time Interval (Min)	Operation Time Interval (Sec)	Cycle Time	Zero Tube Diameter	Zero Tube Length	One Cycle Tube Number	End Time (Min)	End Time (Sec)	
32.1	12.34	10.5	2	30	3	30	1	5	9	26	6	0	

Setting													
Motion Speed	Collect Flowrate	Tube Volume (%)	Start Time (Min)	Start Time (Sec)	Operation Time Interval (Min)	Operation Time Interval (Sec)	Cycle Time	Zero Tube Diameter	Zero Tube Length	One Cycle Tube Number	End Time (Min)	End Time (Sec)	
32.1	12.34	10.5	2	30	3	30	1	5	9	26	6	0	

Run / Stop													
<input type="button" value="Run"/>													<input type="button" value="Close"/>
<input type="button" value="Copy"/> <input type="button" value="Apply"/>													

Fraction													
Status													
Motion Speed	Collect Flowrate	Tube Volume (%)	Start Time (Min)	Start Time (Sec)	Operation Time Interval (Min)	Operation Time Interval (Sec)	Cycle Time	Zero Tube Diameter	Zero Tube Length	One Cycle Tube Number	End Time (Min)	End Time (Sec)	
32.1	12.34	10.5	2	30	3	30	1	5	9	26	6	0	

Setting													
Motion Speed	Collect Flowrate	Tube Volume (%)	Start Time (Min)	Start Time (Sec)	Operation Time Interval (Min)	Operation Time Interval (Sec)	Cycle Time	Zero Tube Diameter	Zero Tube Length	One Cycle Tube Number	End Time (Min)	End Time (Sec)	
32.1	12.34	10.5	2	30	3	30	1	5	9	26	6	0	

Run / Stop													
<input type="button" value="Stop"/>													<input type="button" value="Close"/>
<input type="button" value="Copy"/> <input type="button" value="Apply"/>													

#### 4.2.5 停止餾分: Stop Fraction: Run

Click Stop button

Fraction

Status												
Motion Speed	Collect Flowrate	Tube Volume (%)	Start Time (Min)	Start Time (Sec)	Operation Time Interval (Min)	Operation Time Interval (Sec)	Cycle Time	Zero Tube Diameter	Zero Tube Length	One Cycle Tube Number	End Time (Min)	End Time (Sec)
32.1	12.34	10.5	2	30	3	30	1	5	9	26	6	0

Setting												
Motion Speed	Collect Flowrate	Tube Volume (%)	Start Time (Min)	Start Time (Sec)	Operation Time Interval (Min)	Operation Time Interval (Sec)	Cycle Time	Zero Tube Diameter	Zero Tube Length	One Cycle Tube Number	End Time (Min)	End Time (Sec)
32.1	12.34	10.5	2	30	3	30	1	5	9	26	6	0

Run / Stop												
<input type="button" value="Stop"/> <input type="button" value="Copy"/> <input type="button" value="Apply"/>												
<input type="button" value="Close"/>												

Fraction

Status												
Motion Speed	Collect Flowrate	Tube Volume (%)	Start Time (Min)	Start Time (Sec)	Operation Time Interval (Min)	Operation Time Interval (Sec)	Cycle Time	Zero Tube Diameter	Zero Tube Length	One Cycle Tube Number	End Time (Min)	End Time (Sec)
32.1	12.34	10.5	2	30	3	30	1	5	9	26	6	0

Setting												
Motion Speed	Collect Flowrate	Tube Volume (%)	Start Time (Min)	Start Time (Sec)	Operation Time Interval (Min)	Operation Time Interval (Sec)	Cycle Time	Zero Tube Diameter	Zero Tube Length	One Cycle Tube Number	End Time (Min)	End Time (Sec)
32.1	12.34	10.5	2	30	3	30	1	5	9	26	6	0

Run / Stop												
<input type="button" value="Run"/> <input type="button" value="Copy"/> <input type="button" value="Apply"/>												
<input type="button" value="Close"/>												

#### 4.2.6 關閉餾分: Ctop Fraction:

Click Close button.

Fraction

Status												
Motion Speed	Collect Flowrate	Tube Volume (%)	Start Time (Min)	Start Time (Sec)	Operation Time Interval (Min)	Operation Time Interval (Sec)	Cycle Time	Zero Tube Diameter	Zero Tube Length	One Cycle Tube Number	End Time (Min)	End Time (Sec)
32.1	12.34	10.5	2	30	3	30	1	5	9	26	6	0

Setting												
Motion Speed	Collect Flowrate	Tube Volume (%)	Start Time (Min)	Start Time (Sec)	Operation Time Interval (Min)	Operation Time Interval (Sec)	Cycle Time	Zero Tube Diameter	Zero Tube Length	One Cycle Tube Number	End Time (Min)	End Time (Sec)
32.1	12.34	10.5	2	30	3	30	1	5	9	26	6	0


Run / Stop												
<input type="button" value="Run"/> <input type="button" value="Copy"/> <input type="button" value="Apply"/>												
<input type="button" value="Close"/> <span style="border: 2px solid red; border-radius: 50%; padding: 2px;">Close</span>												

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
	<p>超載, 分離切分使得分離組分成為主成分(富集)後, 再次分離製備.</p> <p>When the fractions were collected, it generally have the following situations:</p> <p>(A) Get a good peak separation: preparative column could be overload to improve efficiency.</p> <p>(B) Being a small component between the two main components:</p> <p>Overload and try to be separated the small component slice · and to enrich it to become a main component, then the enriched preparation is again separated.</p>
(2)	<p>中心切割技術 Central cutting(Splicing) technology</p> <p>分析型分离 → 完全分离组分a的载样极限 → 超载</p>
(3)	<p>邊緣切割技術 Edge cutting (Splicing) technology</p>

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On the Home Page, click Welcome to enter into monitoring interface of the fraction collector. as follows:

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Boot with no display: Check the power input fuse is blown,  
(replacement 10A fuse)

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Ractions connector leakage: Flush fluid valve set with Isopropyl Alcohol solution (valve set may have crystal or dirt accumulated)

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Keypad failure : Contact the manufacturer, non-professionals should not to open the main case (to prevent bodily injury caused

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Solenoid valve problem:

Solenoid valve do not function in fraction when the ON key is pressed, (press the ON button with no action sound) : Contact the manufacturer, non-professionals should not to open the main case (to prevent bodily injury caused by electric shock)

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For users complying with the conditions of storage and use rules, the warranty time is within 24 months from the date of shipment from the manufacturing units. (the seals and vulnerable parts are excluded)

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It is not warrantied for damage from improper operation of equipment, and will have corresponding charges.

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GL-LF096 fraction collector, has been properly packaged at the factory. If any damage of the GL-LF096 fraction collector occurs during transportations, ask statement as soon as possible to the instrument carrier sector for the relevant payment · and notification to the manufacturer for repair or replacement.

### Ordering & Service Information: 訂購和服務信息.

Welcome to our web site to see our currently updated product information and catalog.

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NO.137, JIAZIDONGDAO, DAYI COUNTY, CHENGDU, CHINA 611330

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電 話(Phone) : 028-88201806/88201866, 傳 真(Fax) : 028-88201866,

Web: <http://www.cdgelai.com>, Email: [cdsglp@163.com](mailto:cdsglp@163.com)

The Analytical Based Development Center (ABDC Work Shop)

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Zipcode:407, Add: 5th Floor, NO.641, Fu Shun Road, Shi-Tuen District , Taichung City, Taiwan, R.O.C.

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[http://chromnet.net/Product%20Atlas\\_Preparative%20Chromatography%20Systems.aspx](http://chromnet.net/Product%20Atlas_Preparative%20Chromatography%20Systems.aspx)

Email: [service@chromnet.net](mailto:service@chromnet.net)

LINE:service.abdc, SKYPE: skypeabdc, Twitter: ttabdc, WeChat: wcabdc

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ChengDu GeLai Technology Co.,Ltd.