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Static Axial Compression preparative chromatography Installation & User Manual

靜態軸向壓縮製備層析管柱 Static Axial Compression preparative chromatography column

成都格萊精密儀器有限公司(Chengdu GeLai Technology Co., Ltd)


太陽系分析基礎開發工作室(The Analytical Based Development Center)

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Static Axial Compression technology (SAC)

<ol style="list-style-type: none"> 1. 手動加壓手柄, 2. 行星減速齒輪 3. 減速機機座 4. 防回彈裝置 5. 防回彈可調裝置 6. 層析管柱頭 7. 高分子密封圈 8. 四氟以烯密封圈 9. 多點強制分配器 10. 不鏽鋼過濾器 	<ol style="list-style-type: none"> 11. 層析管柱 12. 層析管柱固定相 8. 四氟以烯密封圈 10. 不鏽鋼過濾器 9. 多點強制分配器 7. 高分子密封圈 6. 層析管柱頭 13. 管柱支撐座 14. 層析管柱支撐腳
<ol style="list-style-type: none"> 1. Manual pressing handle 2. The planetary reduction gear 3. Reducer base 4. Anti-rebound device 5. Anti-rebound adjustable device 6. Chromatography column head 7. Polymer seal 8. Tetra-fluoro-ethylene seals 9. Multi-points distributor 10. Stainless steel filter 	<ol style="list-style-type: none"> 11. The chromatography column 12. The stationary phase 8. TPFE seals 10. The stainless steel filter 9. Multi-points distributor 7. Polymer seals 6. Chromatography column head 13. The column support base 14. The column support legs

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The high density and homogenous distribution of the packing materials within the column bed can be kept in long-term, just necessary for the user to gently tighten manually the re-pressing device on the column top °

It solves the shrinking and collapsing problems during the application courses of the fixed packed column bed, and avoids the cumbersome process of scraping the upper contaminated packing layer.

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SAC and DAC (Dynamic Axial Compression) preparation columns lead the innovation development of preparative column technologies, which can provide maximum sample capacity and lowest cost per unit volume, can provide same or similar separation packing capacity and separation efficiency as that from other high-performance analytical chromatography columns.

D80mm SAC type column, when used in all areas of chemistry are designed with 25cm bed length.

The Gelai SAC column has been widely used in the separation and purification of proteins and peptides and natural products.

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SAC column structure and installation

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SAC column may be used alone, can also be directly connected in series with the same type of column, and can also be pre-equipped with 5 cm long cylindrical Guard column or Online Filter to filter impurities and for the retention of contaminants.

SAC End-device can form high pressure sealing on the inner surface of the column. enabling easy recompressions performed without any tools.

The SAC column system hardware is designed to ensure low initial cost, long life and low cost of updating, but also has a high column efficiency.

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
Compared with other fixed size flange type column, the static axial compression characteristics of Gelai SAC column , eliminates the gap within packed materials, extending the column lifetime.


In the column entrance of the Gelai SAC system, there is a real axial compression mechanism which includes a planetary gear device to ensure to reflect its mechanism advantage in the column packing and compressing process.

Just without tools, gently tightened drive wheel, to push forward the piston in column and realize re-compression and entrance gap elimination function.

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Recompression device of SAC column head

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
Schematic view of the inlet end of the ID80mm SAC column type

The arrow indicate the force and motion direction

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Column support base and support leg

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Each pre-packed SAC column has a high-pressure compressed adsorbent bed that is fully packed with packing material (C18, etc.). There is a movable "column bed driving seat" at the column inlet end, and a fixed one at the column outlet end

The "column bed driving seats" are placed inside column, located near to both ends .

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The direction of column inlet and outlet ends:

- A. Arrows on the label must be the same with the direction of fluid flow.
- B. The wrench flats or grooves of the column outer wall are located near the inlet end.

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Inside the column the SAC axial compression piston can reach sealing against high pressure. The piston can be integrated closely with the column inner wall section between the end surface and "column bed driving seat".

When the column packing is not compressed tightly, the seal ring can slide very easily.

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Manually tighten the drive wheel on the SAC axial compression module will push the axial compression piston toward column bed driving seat of the inlet end。

Manually tighten the drive wheel, to maximize the axial force from piston to the column bed driving seat。

The axial force causes the column bed driving seat of the inlet end to slowly compress packing bed and compensate the packing material reduction to maintain optimal bed packing structure.

SAC 軸向壓縮活塞和管柱床傳動座中均有液流分配器，用以把溶劑均勻地引入填料床的入口端和引出出口端。這些通道經過了最優化，以獲得最小的帶狀鋪展。HPLC 系統的 1/8" 管經由帶紋軸襯和金屬套環與軸向壓縮活塞連在一起。Both the SAC axial compression piston and the column bed driving seat have fluid distributor that are used to uniformly introduce the solvent into the inlet end and out of outlet end of the packed bed. These passages have been mostly optimized to obtain minimal strip spreading.

The 1/8 " tube of the HPLC system are connected together with axial compression piston via patterned bushing and metal ferrule.

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SAC high pressure seal is a passive seal, similar to the HPLC pump piston seal. The seal is a concave casting ring with oblique spring, and with one open end and the other closed end, and the open end always towards to the inside of the column .

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When the column packing is not pressured, the oblique spring ring within the seal still has enough tension to keep the seal in contact with the column wall at low pressure status . and also without affecting the seal to slip when manually tighten the drive wheel .

When the column is pressured, the solvent enter the ring from the open end , and as the pressure increases, it press the seal to closely contact with the column wall and the piston body. In the HPLC system pressure, external forces can also play a role in sealing.

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Warning: Effective sealing of SAC column depends only on the passive sealing mechanism. It only needs by hand to fasten the end device with no need to use greater force.

Fastening with tools will not enhance sealing, but on the reverse, will

damage the column and packed bed. This is particularly important for ID80mm type SAC system.

In the compressing packed bed, a planetary gear unit can provide enough mechanical advantage.

*** Prohibit the use of tools, only hand tight wheels ***

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The easiest SAC column configuration includes a Online filter and column end devices.

Installation SAC Column:

The column was removed from the box. The column should be covered with plastic covers to protect threads on the column during transportations.

(Note: The plastic cover is not included in pre-packed column).

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Take the end devices from the packaged box, and check the following parts is available:

An axial compression module (Note: The packaging covers are not included in pre - packed column)

An outlet axial compression nut / column seat and support(including four supporting legs)

An inlet axial compression piston(with 4 hexagon socket screws, corresponding to 1/8" inlet side)

An spare inlet axial compression piston

(Note: The packaging covers are not included in pre - packed column)

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Each box is also provided with parts for connection between the column and the HPLC system.

These connection parts include a female - female union for 1/8 " outer diameter (OD) pipe and a reusable plastic two-way type device connector, a 50cm long × 1/8" OD × 0.062 " PEEK tubing with a 1/8 "nut and ferrules on the ends.

One end of the tube is inserted into the union, the other end is connected to column inlet.

With the union, the reusable plastic device connects the column and HPLC system together.

There are also two 1/8 "nuts and two ferrules spare parts.

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Check each axial compression piston, ensure that seals included, and the open side of seal (you can see the spring) faces against the sealing drive seat (Note: The above parts is not visible in the pre - packed

column).

Slide the axial compression piston from and into the column outlet.

The arrow on the column label indicate the direction of the liquid path, with the outlet axial compression piston mounted on one end farther away from the entrance of the column.

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Push with more strength the axial compression piston into the column · then mount the seal ring into column °

Install outlet compression nut and column support base, and put the four support legs into the corresponding nut and tighten the screw to complete and form a four-base support, enabling the stand status of the separation column.

Manually screw to tighten the outlet axial compression nut to column outlet, and

Rotate manually to tighten the retaining screws to ensure a solid foundation.

The column should be placed in a horizontal plane.

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Check the column inlet wall, make sure there is no residual particule material, which is particularly important to prevent abrasion of seal ring. Slide the inlet axial compression piston into column inlet. With

more strength by hand to push the axial compression piston into column, then mount the seal ring into the column.

Connect the inlet axial compression piston to female - female union with pipe.

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First, counterclockwise to unscrew and remove the drive wheel of axial compression module. Along the threads to tight the union. Pass the connecting pipe through the center hole of the axial compression module for connection.

Slowly by hand to tightened axial compression module to the column inlet, until it feels tighted (no need to use a lot of strength).

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Prohibit the use of tools !

Man power can ensure appropriate axial force, and to prevent too tight that will damage the column structure.

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Slowly and clockwise to tighten the drive wheel, until it feels its embedded piston has contact with the inlet of axial compression

piston. While touched, there will be a slight, not continuous, and perceptible increase of resistance, that will resist the rotation of the drive wheel, and then tighten the wheels once circle again.

Connect the assembled column to HPLC system with the provided nuts and ferrules. Confirm that the fluid direction is in consistence with that indicated by the arrow on the column label. See NO.38-3060,38- 3061,1533, U-175 or U-176 catalog items for All connection parts .

A plastic and reusable connector (see catalog item NO. FTF-110) is recommended for connecting the column outlet to the system.

On the first time to use SAC column, please do the following steps in for mobile phase equilibrium.

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Before each sample loading, the column mobile and stationary phases must be equilibrated.

Ensure that the mobile phases through the column are compatible with column packing material · the stationary phase.

See Chemistry Guide for reference.

According to the mobile phase isocratic or gradient program · pump sufficient volume of mobile phase through column until the phases reach equilibrium .

Use detector to monitor to see if the phase equilibrium reached .

The mobile phase volume needed depends on the phases used .

For generic cases , 4 times column volume is necessary .

The packing volume of the ID80mm * 25cm (length) SAC column is about 1250ml.

When baseline drift reaches the minimum , and the retention volume can be repeated in successive test runs, the column has reach the equilibrium state.

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Counter-clockwise rotate to remove the axial compression nut from both ends of the column.

Before the axial compression module is removed , Counter-clockwise to rotate the drive wheel , until its embedded piston is no longer in touch with the inlet axial compression piston , and then unscrew to remove the axial compression module.

Draw the axial compression piston from both ends of the column.

In case the column has been separated from the system, then the embedding of a 10-32 or 1/8 "bushing into the axial compression piston will help to easily pull out the piston from the column by holding the bushing .

Take the new packaging cover on the column down, and follow the steps in part 3 to install.

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If you want to get down for preservation of the SAC column from HPLC system, follow these steps:

When the column is still connected to the system, equilibrate with the solvent for preservation. (See part 4, process according to the column equilibrium instructions in order to maintain stability. See Chemistry Guide for reference)

The column was then removed from the system (see part 5) and the package cap covered on both ends of column.

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Before the material enters the SAC column, the Online filter or Guard column, can be used to filter impurities and contaminants.

Also, since the preparative HPLC separation conditions must be

Since the preparative HPLC feed must be well pretreated to maximize yield with minimized solvent and media cost, the use the online filters prefers.

The online filter membranes can be replaced easily.

The online filter has very small volume factor and tiny fluid diffusion

effect that is particularly suitable for the protection in front of the preparative HPLC columns.

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The guard column is filled with the same packing material as that of separation column. It can filter impurities and retain pollutants, but the guard column has definite volume, that exhibits fluid diffusion effect to certain extent.

And hence · the feed (raw material) should be well prepared as possible in advance · and only adopt guard column in cases that the feed(raw material) preparation is complex and has high uncertainty of impurities and contaminants 。

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8.1 Installation of column with guard column :

To installation a column with guard column, a packed column, a packed guard column, and column end devices box is necessary 。

The guard column is filled with the same packing material as that of separation column.

For adding a guard column to existing column, please refer to part

9.

Remove the packaging covers from column and guard column .

Take the end devices box from packaging case.

Check the following components are available:

- ✓ 軸向壓縮模組件一件
- ✓ 出口軸向壓縮螺帽一件/管柱支架 (帶 4 條腳)
- ✓ 入口軸向壓縮活塞一件(4 個內六角螺釘, 1-32-1/8 進口端相對應)
- ✓ 出口軸向壓縮活塞一件(2 個內六角螺釘, /4-28 出口端鄰近)
- ✓ 軸向耦合活塞一件(活塞有兩個密封圈)
- ✓ 前置保護管管柱適配器插套一件
- ✓ 備用入口軸向壓縮活塞體一件
- ✓ an axial compression module
- ✓ a outlet axial compression nut and column base support(with four legs)
- ✓ an inlet axial compression piston (4 hexagon socket screws, 1-32-1/8 corresponding to inlet end)
- ✓ an outlet axial compression piston (2 hexagon socket screws, /4-28 near the outlet end)
- ✓ an axial coupling piston (piston with two seal rings)
- ✓ a guard column adapter plug
- ✓ an inlet axial compression piston body

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Each device box contains all the parts(the end devices box description in part3) required to connect a column with guard column to HPLC system

Check each piston . confirm the correct seal ring direction is located

(see Part 3),

Also check the column and guard column wall to remove possible presence of micro-slag .

Slide the outlet axial compression piston into the column outlet (see Section 3).

Install outlet axial compression nut and column bracket (see section 3).

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Slide the coupling piston into the column inlet, (see Part 3, similar to the inlet axial compression piston installation).

Set guard column adapter plug on the guard column and carefully fix it firmly on the column mouth.

Do not over force. Ensure that the thread on the adapter plug can be screwed more easily. If screwing blocking encountered, screw it down, and try again after cleaned. Do not force and prohibit using tools!

Warning: The adapter plug is very important. Note to protect it against falling.

Use pipe to connect inlet axial compression piston and unions.

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Lead the union and the connecting pipe through the central hole in the axial compression module, and hand-tighten the screw to connect the guard column adapter plug. Be careful not to be too forced.

Slowly and clockwise to tighten the drive wheel, until it feels its embedded piston has contact with the inlet of axial compression piston (see part 3), then tighten the drive wheel one circle.

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To standing for a few minutes is necessary. It can purge out air between guard column and column, to avoid the formation of dead volume that may result in column performance degradation.

Tighten the drive wheel one circle turn.

Connect the SAC column with guard column to HPLC system, as in part 3.

In accordance with the described steps in part 4 to equilibrate mobile and stationary phases for the column before sample loaded.

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8.2 Replacing guard column：

Guard column replacement cycle: According to the complexity of sample preparation for the actual impurities and contaminants situations · to determine the required replacement frequency · Periodically and/or conditionally replace the guard column to protect against chemical contamination.

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Replacing Procedure:

Release HPLC system to zero pressure · then

remove the axial compression module, as described in Part 5.

Pull out the inlet axial compression piston from the guard column.

If the column has been separated from the HPLC system, then a 1. / 16 "-10-32 or 1/8" bushing could be embedded in the piston, to let easily grasping the bushing by hand to remove the piston.

Remove the old guard column, replace it with a new one. Confirm the direction of the arrow of the guard column label is consistent with the direction of the fluid.

Slide the inlet axial compression piston into the new guard column, and again to conduct the connection, exhausting air entered and other steps, as described in Section 8.1.

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8.3 Prepend prepend

With a "guard column-column" coupling devices, a guard column can be added in front to a pre-packed column. The diameter and packing material of the guard column must be the same with that of the connected column.

The "guard column-column" coupling devices include a column adapter of guard column and an axially coupling piston.

Check if all the parts are available, and check if the seal ring on the axially coupling piston is placed on correct direction.

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Adding the guard column to column:

Separate the column inlet from the HPLC system.

Counterclockwise rotate the drive wheel of the axial compression module by hand · until that the piston within the axial compression module is no longer in contact with the inlet axial compression piston. Then, counter-clockwise rotate to take the axial compression modules out.

Pulled out the inlet axial compression piston from the column.

Slide the axial coupling piston into the column inlet. That is step to substitute the removed axial compression piston.

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Set the guard column outlet to the axial coupling piston ·

Confirm the direction of the arrow of the guard column label is

consistent with the direction of the fluid.

Slide the inlet axial compression piston that is take from column, into the inlet of guard column.

Set the "guard column adapter plug " on the guard column, hand-tighten onto the chromatography column.

Connect the inlet axial compression piston to union with tube °

Connect axially compressive module, as described in part 8.1.

Exhaust entered air, as described in Section 8.1.

Reconnect the column inlet to HPLC system.

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More than one column could be serially connected with a minimum length of stainless steel tubing (0.062 "ID type).

Connect the outlet of one column to the inlet of another column with a suitable 1/8 "hand - tightened joint (or plastic reusable connector).

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When a new column installed, or when not in use for a long time, it is usually required to evaluate its performance using standard test methods. See the specific test conditions detailed in the Handbook of Chemistry.

For precaution and usage of SAC column · please refer to chromatography column related Handbook of Chemistry

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The Interspace in the packing may lead to reduced performance of the column.

In order to ensure long life of the column, it is required to periodically to release the HPLC system pressure to zero, and then manually rotate to tightened all the axial compression parts of the SAC column. That can maintain optimal packed bed layer structure, and minimize formation of packing bed Interspace.

Do not tighten SAC axial compression device when the column of the HPLC system is under pressured status, since the seal rings can not slide easily.

So, the HPLC system pressure must be completely released before tightening.

The SAC column system is designed to ensure that hand - tightening is sufficient.

Excessive tightening by tool may damage the crucial internal parts.

Do not use tools to tighten SAC axial compression device.

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Importantly:

If leaking around seal rings of the SAC column occurs, the seal rings should be replaced, not resolved by increasing the axial pressure. That is, it is invalid by rotate to tighten the drive wheel of axial compression module.

Keep connecting pipe to the HPLC system to minimum length to avoid and eliminate dead volume that may be left in the device and connectors.

When connecting to system, do not over tighten the bushing and the metal ferrule.

Excessive tightening by tool may damage the axial compression piston.

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When connecting guard to the column, be sure to follow the instructions, to purge the introduced air out of the axial coupling piston, or it may produce dead volumes and column performance degradation.

Use only HPLC- grade solvents, which need to use helium degassing or vacuum filtration.

When all of the mobile phase was filtered, the Online—Filter with 0.22um membrane is recommended to filter micro-particles on line (using membrane material tolerable to HPLC solvents).

It is a fundamental demand for solvents with dissolved salts or for

solvents that may support bacteria multiplication.

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Adjust the column to equilibrium state before the sample loading at the beginning of chromatographic separation,

Do not exceed the pressure limit of the column.

The pressure limit of the pre-packed ID80mm type SAC column is about 2000 psi (140Bar). The pressure limit of the ID80mm type SAC column hardware is 4000 psi (280Bar) . Once the overvoltage, the first inlet seal rings is usually the first part to be damaged.

The Pulse damper reduces vibration. Vibration is the main reason of causing reduced column life. It is especially important when the solvent flow very fast.

The column performance deterioration due to chemical contamination can usually be recovered with a cleaning method, but the solvents have more concentration than that used in routine chromatography.

For more information, please refer Handbook of Chemistry.

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12.1 Replacing the Online Filter

Completely release the column pressure to zero.

Remove the axial compression module.

From the inlet axial compression piston with online filter, sequentially removed the four hexagon socket screws, the seals and seal retainer.

Install a new solvent online filter to the new axial compression piston, then according right direction to install seal ring and seal ring holder, and tighten the hexagon socket screws.

Place the inlet axial compression piston to the inlet of column, and re-install the axial compression module.

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12.2 Replace axial compression piston seal

Completely release the column pressure to zero.

If the inlet device leaks · the axial compression module have to be removed.

If the outlet device leaks · the axial compression module nut / column support base(the base cap)have to be removed as following:

Separate the column outlet apart from the HPLC system.

Counter-clockwise rotate thumb screw bar to remove it, while the column is supported by hand. In counterclockwise direction, remove the axial compression screw cap(the base cap) by hand. Do not need to remove the nuts on the cap · for it forms a solid base of four.

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Remove the axial compression piston from the column.

Unlock and remove the hexagon socket screws (two on outlet piston, four on inlet · these bolts link the seal ring retainer and piston body together) · to separate the seal ring retainer and piston body apart ·

Remove the old and replace the new seal ring. The open end of the seal (with visible spring) should face to the seal ring retainer.

Replace seal ring retainer and reinstall the device.

The replacement of the seal ring of the axial coupling piston is similar.

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備用裝置和連接管 Backups and the connection pipe

配件編號	描述 Part Number / Description
38-3060	PEEK, 1/8 OD x 0.062"OD 型 x 50cm PEEK, 1/8 OD x 0.062 "OD x 50cm type
38-3061	PEEK, 1/8" OD x 0.062"OD 型 x 50cm PEEK, 1/8 "OD x 0.062" OD x 50cm type
01-0140	手動旋緊接頭, 帶 1/8"-螺紋, 配 1/8"x 0.062"OD 連接管, 2 只。 Manually tighten fitting with 1/8 " - thread, with 1/8" x 0.062 "OD connection tube , two pieces.
01-0142	不銹鋼金屬套環, 1/8", 10 只。 Stainless steel ferrule, 1/8 " , 10 pieces
01-0150	不銹鋼兩通, 配 1/8"OD , 1 只。 Stainless steel two - way with 1/8 "OD, 1 pieces.
FTF-110	塑膠可重複使用連接件, 配 1/8"連接管.與不銹鋼型裝置用螺紋相配

	合 Reusable plastic fittings, with 1/8 "connection tube to match with stainless steel screw-type device.
--	--

備用活塞配置 (配密封圈和固定器)
Spare piston configuration (with seal rings and retainer)

配件編號	Part Number / Description
7002-211	80mm ID 型系統的入口軸向壓縮活塞裝置 Inlet axial compression piston device for 80mm ID System
7002-059M	80mm ID 型系統的出口軸向壓縮活塞裝置 Outlet axial compression piston device for 80mm ID System
7002-060M	80mm ID 型系統的軸向耦合活塞裝置 Axial coupling piston device of for 80mm ID System

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For more application specific information and data for column and guard column when ordering, see Handbook of Chemistry.

Ordering & Service Information: 訂購和服務信息.

Welcome to our web site to see our currently updated product information and catalog.

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Static Axial Compression preparative chromatography Installation & User Manual

NO.137, JIAZIDONGDAO, DAYI COUNTY, CHENGDU, CHINA 611330

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電話(Phone) : 028-88201806/88201866, 傳真(Fax) : 028-88201866,

Web: <http://www.cdgelai.com>, Email: cdsglp@163.com

The Analytical Based Development Center (ABDC Work Shop)

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Zipcode:407, Add: 5th Floor, NO.641, Fu Shun Road, Shi-Tuen District , Taichung City,
Taiwan, R.O.C.

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http://chromnet.net/Product%20Atlas_Preparative%20Chromatography%20Systems.aspx

Email: service@chromnet.net

LINE:service.abdc, SKYPE: skypeabdc, Twitter: ttabdc, WeChat: wcabdc

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ChengDu GeLai Technology Co.,Ltd.