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UV-3292 UV spectrophotometric detector operating manual

成都格萊精密儀器有限公司(Chengdu GeLai Technology Co., Ltd)
太陽系分析基礎開發工作室(The Analytical Based Development Center)

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The UV-3292 UV spectrophotometric detector is implemented with full digital output design。The analog detector signal is converted with zero adjustment by the microcontroller and 24-bit AD converter, with the processed results output to the RS232 interface.

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The grating of the detector is fixed directly to the motor shaft of the stepping motor。The sinusoidal conversion by the microcontroller control the motor to wavelength selection that replaces the traditional sine bar mechanism, so that the reliability and stability of the system has been further improved.

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Wavelength calibration is built in system boot steps.

Programmed wavelength control.

Stopped-flow spectral scanning capability

The flow cell is design with Z-shaped and thermostatic features, that improving the instrument noise and drift minimization.

Double beam optical design, further reducing the environmental effect on signal drift.

The light path can be adjusted by changing the optical path of the flow cell and the corresponding system parameters. It can be easily converted for analytical, semi-preparative, large preparative HPLC.

Analog output port for chromatogram acquisition

System controlled by chromatography workstation through RS232 / RS485 / USB interfaces.

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Adopt Deuterium lamp (D₂ lamp) for light source, plane grating for spectral resolution element, and photodiode (photocell) for a signal-receiving element. The next picture shows the optical path diagram of the UV spectrophotometric detector:

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Figure 1 Optical path diagram of the UV spectrophotometric detector

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1.Deuterium lamp, 2.Slit, 3.Filters, 4.Concave mirror, 5.Plane, 6.grating,
7.monochromator, 8.Spectroscope, 9.Flow cell,
10.Photodiode (for sample slot), 11.Photodiode (for reference slot)

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A brief description of the Optical path diagram:

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As shown in Figure 1, the deuterium lamp(1) light beam emits from the slit(2), through filter(3), then focused by the concave mirror(4) to the plane grating monochromator(5), thereby obtaining a light beam of a specific wavelength. This beam is divided into two portions by a dichroic mirror(6). One portion of the beam is transmitted to a reference photodiode(9) to generate a reference signal, and the other is conducted to pass through the flow cell(7), that exposure to the photodiode (8), generating a measuring signal.

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The front panel of the UV spectrophotometric detector if composed of three parts:

Information display area, Keypad area and the optical unit components. As in Figure 2, The Display and Keypad sections are used for displaying and controlling detector status and operation parameters。

The UV detector pipe connection, flow cell installation and maintenances, can be reached directly from the front panel optical unit.

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1) Information display area 2) Keypad 3) Optical Unit Assembly

3.3 流通池 Flow Cell

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Figure 3 Flow Cell areas and Flow cell

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As shown in FIG. 3, 4), Liquid inlet is connected to chromatography column. The pipe must be inserted as close as possible to the flow cell to prevent from producing dead volume.

As in FIG. 3, 5), Liquid outlet is connected to a back pressure

device, then goto the waste or fraction collection device.

As in FIG. 3, 1), The flow cell can be removed for cleaning and maintenance by loosen the two fixing screws of the flow cell.

See the instructions in "5.1 Cleaning flow cell" and 5.2 "Removing the flow cell for cleansing"

Be sure to install the flow cell with the fixing screws fixed firmly to prevent from leakage of light.

3.4 顯示區和薄膜鍵盤 Membrane keypad and display area

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Figure 4 Display area and Membrane keypad

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By operating the up and down paging keys, the Information display area displays pages for the current settings for the wavelength, signal output value, D₂ lamp energy and accumulated operating time, serial number and operating

functions menu, logo menu.

Membrane keypad is used input data, erase the data, and to achieve detector auto-zero and program operation control and other operations .


Program running of the detector can also be controlled by the computer workstation software .

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The LCD screen displays the main interface following the Instrument start up, as the display status:

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In this case the detector wavelength is 240nm. ">" means that the wavelength can be set with digital and Enter keys;

-0.0000AU is the signal output value (AU), and SN: The instrument manufactory serial number;

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Through the digital keys and Enter key to set the displays the time constant of the filter: RESPONSE 1.00 (seconds) · that is usually set to 1 (sec). Press Enter and then to set analog signal range :AUX RANG 1.00 (AU / V) · that usually is set to 1.00 (AU / V)

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Figure 5 UV spectrophotometric detector rear view

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As shown in Figure 5， 1) RS-232 communication interface (9-pin connector). Through the RS-232 cable connected to PC, communication, the exchange of data with the PC is completed. Also, a pre-installed PC workstation software can control the UV spectrophotometric detector.

3) The DC24V host shell fan is used to dissipate and balance heat form the host shell transformers and detector deuterium lamp.

2) 220V power input, built two fuses (250V1A), one of which is spare for convenience of replacement。

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如使用模擬口，還需通過 D/A 轉換器 Uv-3260 轉換後再連工作站轉接器後方可與連接電腦。
For UV-3292/UV-3293 to communicate with PC, the computer must have a free serial RS232 port or a extended serial interface, such as the 『USB to RS232』 convertor. Find the RS232 cable to the accessory kit, and connect firmly one end of the communication cable to the rear panel RS232 communication interface, and the other end to PC serial port(or USB to RS232 convertor), and then correctly set up the corresponding port position on the workstation configuration.

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1. Use the instrument parts and accessory to connect detector pipes, paying attention to correctly locate the flow cell liquid inlet in lower position and outlet higher position that may avoid the accumulation of bubble inside the cell.

2. Insert the end of power cable into power input socket of the instrument rear panel, and plug the other end into the power source socket, then turn ON the power switch on the right side of the instrument rear panel, 0 is Off, 1 is On.

3. Self-diagnosis and wavelength calibration is built in system boot steps.

4. After successful Self-diagnosis, the instrument is in standby mode, with D2 indicator lighting (Deuterium lamp is lighted), then you can set the Wavelength (nm) with the keypad, the Filter time constant RESPONSE:1 (usually set to 1 second), and analog signal range AUX RANG: 1.00 (usually set 1.00 AU/V).

5. After the instrument have warmed up to stable readings, press Auto Zero key to zero (minimize) the detector baseline reading.

6. Start to detect baseline and chromatograms.

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Special statement : The opening of this product casing by non-professional or authorized personnels is prohibited .

It will cause invalid Warranty conditions.

5.1 Cleaning the flow cell Cleaning the flow cell

Greater baseline noise and low sensitivity may be caused by contaminated flow cell. Rinse and flush the flow cell according to the following procedures, and in most cases it is effective.

a) Use one of the following solvent to wash the flow cell:

Sodium dodecyl sulfate, 1m HCL, 1m NaOH, Ethanol, Acetone.

b) Injected the solvent flow cell by a syringe, stay the solvent inside for about 5 minutes.

c) Rinse thoroughly with water, then use with a gentle stream of pure Nitrogen air to rush to dry.

5.2 拆卸流通池清洗 Removing the flow cell for cleansing

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Figure 11. A cross-sectional view of the flow cell

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a) Use the 3mm hex wrench in the flow cell spare parts to loosen the outer nut.

b) Removed the lens-equipped PEEK holder with a clamp, orgently knock it out onto a flat, clean and soft desk surface.

The Lens is embedded in the PEEK holder with PTFE sealing ring, demanding that each disassembly and assembly of the flow cell should replaces the seals.

c) Remove the lens, wipe with a soft cloth or select the appropriate solvent for ultrasonic cleaning. Be careful not to touch the clean lens with your fingers.

d) Reverse the above procedure in order to Reinstall the flow cell, and make sure not to block the optical path by the PTFE seals.

e) Tighten gently and carefully the outer nut with a wrench to protect

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UV-3292 / UV-3293 UV spectrophotometric detector used Deuterium lamp as light source, and its warranty life is 1000 hours.

The actual service time of the Deuterium lamp depends on various factors such as the number of ignitions, average lighting time and noise and sensitivity requirements.

Check the light energy, both the sample and reference light intensity values .

The reference value being the light intensity value measured from the reference channel, and it can be used to check the quality of the lamp.

Periodic inspection of reference values (virtual flow cell, a wavelength of 240nm) is necessary , in particular, check this value when the detector appears a larger noise increasing or sensitivity decreasing operating states.

When the reference value is close to or less than 0.1 or Deuterium lamp lighting time exceeds 1000h, it must be replaced with new deuterium lamp.

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For users complying with the conditions of storage and use rules, and if the products and parts have manufacturing quality defect and does not work correctly, the manufacturer should be free of charge for the user to repair or replace.

The warranty time is within twelve months from the date of shipment from the manufacturing units.

For consumables, such as the Deuterium lamp, the time of warranty for is accumulative. (The Deuterium lamp, the warranty time is 2000 hours)

Also, it is not warranted for damage from improper operation of equipment.

Ordering & Service Information: 訂購和服務信息.

Welcome to our web site to see our currently updated product information and catalog.

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NO.137, JIAZIDONGDAO, DAYI COUNTY, CHENGDU, CHINA 611330

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電話(Phone) : 028-88201806/88201866, 傳真(Fax) : 028-88201866,

Web: <http://www.cdgelai.com>, Email: cdsglp@163.com

The Analytical Based Development Center (ABDC Work Shop)

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Zipcode:407, Add: 5th Floor, NO.641, Fu Shun Road, Shi-Tuen District , Taichung City, Taiwan, R.O.C.

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http://chromnet.net/Product%20Atlas_Preparative%20Chromatography%20Systems.aspx

Email: service@chromnet.net

LINE:service.abdc, SKYPE: skypeabdc, Twitter: ttabdc, WeChat: wcabdc

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ChengDu GeLai Technology Co.,Ltd.